



## The first profitable silver recovery device designed for the smaller darkroom...that really works!

For years smaller darkroom operators have been pouring valuable silver down the drain with their used hypo/fix solution. Mostly because they don't know this precious silver can be recovered; how valuable it really is; or they consider recovery too complicated and expensive.

The Silver Magnet (SMT-20) is a unique, patented electronic silver recovery module that can recover up to a troy ounce of silver per week from your used solution, simply and economically. When full, it holds 30+ troy ounces of 90 to 95% pure silver.

No one is too small. You can collect up to \$150.00 at \$5.00 per troy oz. from less than 100 gallons. And the Silver Magnet is priced so low, it can pay for itself in just a month of use.

So why throw money away? With the continually rising prices of film and photographic paper, you just can't afford to anymore.

**START PROFITING TODAY.  
START USING THE SILVER MAGNET!**

[silvermagnet@nanianphoto.com](mailto:silvermagnet@nanianphoto.com)

# Easy as 1,2,3

# Silver Magnet™

## SMT-20 used world wide by Photographers, Dentists, Printers, Veterinarians, and Chiropractors

**1** Connect power supply to pick-up cell. Plug in power supply and submerge cell in hypo/fixer bath.

**2** When cell becomes full (30+ ounces) disconnect from power supply and send to refinery.

**3** A simple process you can repeat over and over. Watch your profits grow and darkroom costs go down.

The complete system consists of only 2 parts. A reusable power supply that plugs into any 110-120VAC outlet (220VAC model available); and a replaceable pick-up cell.

The pick-up cell is completely transparent, so you can actually see your silver being collected without removing it from the tray or batch tank containing the used solution.

The cell can be left in your solution during processing. When used in this way, it can actually extend the useful life of your hypo/fix, giving you an additional cost cutting benefit.

When full, just put pick-up cell portion into provided leak-proof envelope, addressed to our refinery. You may request to purchase a new cell at this time also. You will be sent a check for your silver, at market value, and your new cell, if requested.

A patented product with years of use in the field by amateurs and professionals alike, the Silver Magnet (SMT-20) has been proven to be an extremely profitable method of silver recovery for the smaller darkroom.

### **SPECIFICATIONS:**

**Recovery Method:** Electrolytic

**Pick-Up Cell:** 7-3/4" x 3" x 1-1/2"

**Power Module Size:** 2" x 2" x 1-1/2"

**Indicators:** Light Emitting Diodes

**Recovery Capacity:** Approx. 3 lbs. Silver

**Recovery Rate:** Up to 1 troy ounce per week

**Power Input:** 115VAC 60Hz 5 watt, 220VAC 50Hz 7 watt

Euro Version

**Power Output:** 2 VDC 30MA

**Connection:** Chemical Resistant

Shielded 3.5 mm plug and line jack

**Catalog # SMT-20**

